



External Protocol for submitting DNA for Illumina sequencing genomes

Version: 1

[Summary](#)

[Reagents and Materials](#)

[Protocol](#)

Summary:

This protocol is for the submission of DNA for the generation of paired-end libraries for sequencing on the Illumina MiSeq platform. DNA is quantified using a fluorometric-based method and diluted to 0.2 ng/ μ L. Libraries are prepared using the [NexteraXT Library prep kit](#), which uses an enzymatic reaction called tagmentation to fragment the DNA and add adapter sequences. After libraries are prepared they are quantified with another low-diversity library and run together on the MiSeq. The “Preparing Libraries for Sequencing on the MiSeq” (part 15039740, Rev. D) protocol was used to prepare libraries with a final load concentration of 5.5 pM, spiked with 85% low diversity library, 13% genomic library and 2% PhiX. FASTQ files are distributed to the client when the 2 x 250 bp sequencing completes.

Reagents and Materials:

Reagent/Material	Vendor	Stock Number
Aluminum foil, sterile	VWR	89049-034
Full-skirted PCR plate (twin.tec)	Eppendorf	951020401

Protocol:

1. In a clean hood, add at least 15 μ L of DNA to each well of the twin.tec plates. Seal with sterile foil seal. Clearly label plates with PI, Reference ID, and date.

IMPORTANT: Seal plates very well to reduce evaporation and cross contamination between wells.

2. Use these [shipping directions](#) to prepare the PCR plates for shipping.

3. Fill out submission [form](#).
4. Send a list of the genome names and shipment tracking information to msmblcore@umich.edu.
5. Ship on dry ice or with ice packs to:

Host Microbiome Initiative
University of Michigan Medical School
Internal Medicine/Infectious Diseases
1500 MSRB1
1150 W. Medical Center Drive
Ann Arbor, MI 48109-5666

Please include reference ID on package documentation.